



AJ 5807 NEONATAL SYRINGE FEEDING PUMP



A sterilized, DEHP-free and latex-free neonatal feeding solution

Using our syringe feeding pump is completely safe, no need to worry about the DEHP using the DEHP free syringes

Simple & affordable, with optional multi-angle pole clamp

Simple & safe operation, use optional multi-angle pole clamp to ensure nutrition, e.g. (fat in the milk) goes to baby first

Convenient feeding time setting operations

Easy and quick feeding time settings (10 / 20 / 30 / 40 / 50min)

Automatic calibrate to different syringe brands/sizes



Use optional multi-angle pole clamp to ensure nutrition, e.g. (fat in the milk) goes to baby first goes to baby first

- 1 Ergonomic touch panel keypad for easy operations
- 2 LED indicator for clear view of the default infusion time settings

Works perfectly with our compatible fluid warmer FW-300 to feed nutritional to neonatal at the right temperature

- 3 LED indicator for a clear view of the power source used
- Audible & visual alarms
- 4 "NEAR END", "COMPLETE" "OCCLUSION", "LOW BAT"

Specifications

Power supply	AC 110V / 60 Hz or AC 220V / 50 Hz
Motor type	Stepper motor driven linear actuator, pulsed motion for easier digest
Applicable syringes	20 ml, 30 ml, 50 ml, 60ml, 100 ml disposable syringes; Compatible with major syringe brands
Time settings	10 / 20 / 30 / 40 / 50min
Infusion accuracy	± 20 sec
Flow rate	0.4 ml/min (20ml at 50minutes) to 10 ml/min (100ml at 10minutes)
Audible & Visual Alarms	<ul style="list-style-type: none"> NEAR END COMPLETE OCCLUSION LOW BAT
Battery	Built-in rechargeable Ni-MH battery. Min. 4 hours of continuous operation at 60 ml/h with fully charged battery
Included accessories	Standard power cable, pole clamp is optional
Dimensions	115 (D) x 366 (W) x 120 (H) mm / 4.53" (L) x 14.41" (W) x 4.72" (H)
Net weight	1.8 kg / 3.96 lbs