

## Product Specifications

Infusion Rate Range	0.1 - 1200 ml/h (or up to 200 drops/min)	
Infusion Modes	ml/h drop/min	volume/time* patient/weight*
Volume To Be Infused	1 - 9999 ml	
Bolus Rate Range	1 -1200 ml/h	
Bolus Volume	1 - 9999 ml	
Accuracy	± 2% on the pump only (± 5% with I.V. set)	
Occlusion Pressure Levels	High, Medium, Low, Lowest 0.04 MP a -0.16 MP a / 5.8 PSI ~23.2PSI	
KVO (Keep Vein Open)	0.1 ml/h - 3 ml/h (depends on the flow rate)	
Built-in Infusion Modes Memory	10 preset infusion modes memory	
Built-in Multi-languages Display	English, Spanish, French, Italian (Optional:Turkish, Russian)	
Alarm Volume Adjustment	High, Medium, Low	
Alarms	DOOR OPEN, OCCLUSION (up & down), POWER SOURCE CHANGE, START REMINDER, INFUSION COMPLETED, I.V. SET DISENGAGED, EMPTY BAG, INITIALIZATION ERROR, DRIVE MECHANISM ERROR, AIR IN LINE, BATTERY LOW, BATTERY DEPLETED.	
Event Log	400 + records of settings, running status, alarms with date & time stamps	
Safety Protection	Anti-free flow protection, Audible & Visual Alarms, Automatic keypad lock, Password protection	
Power Supply	AC 100 - 240V 50/60 Hz; DC 12V	
Battery	Built-in rechargeable sealed NiMH Pump will operate 3 hours min. at 25 ml/h with fully charged battery	
Safety Approvals	CE 0197, IEC 60601 - 1 - 8: 2003 type CF, IPX1	
Dimensions	105(W) x 172(D) x 240(H) mm / 4"(W) x 6"(D) x 6"(H)	
Net Weight	2.9 kg / 6.4 lbs	

\* This feature only available on the AJ5808A.



**Angel Canada**  
 Angel Canada Enterprises Ltd.  
 301-2999 Underhill Ave. Burnaby, BC Canada. V5A 3C2  
 Tel: +1(604) 568.0936 / +1(604) 568.6977 Fax: +1(604) 677.5295

[info@angelpump.com](mailto:info@angelpump.com) / [www.angelpump.com](http://www.angelpump.com)

**Angel China**  
 Shanghai Angel Electronic Equipment Co., Ltd.  
 4/F Building 5, No. 1135, Baotou Road, Shanghai, China. 200438  
 Tel: +86(21) 6532.1221 Fax: +86(21) 6506.7930



# AJ 5808 / AJ 5808A Volumetric Infusion Pump



ID:  
32268







### AJ 5808 INFUSION PUMP WITH SIMPLE KEYPAD

AJ5808 is simply to use, easy to program. It's an economic solution for reliable and safety infusion in any hospital and care centers.



### AJ 5808A INFUSION PUMP WITH NUMERIC KEYPAD

AJ5808A can be programmed in four different infusion modes. It comes with drug dose calculator enabling users to setup the infusion therapy based on patient weight and drug dose safely and quickly.



• Easy to program numeric keypad.



• Alarm & warning beacon.

### MAIN FEATURES

- Micro / Macro infusion rate setting from 0.1 to 1200 ml/hr.
- OLED screen provides a more vivid reading than traditional LCD screens even from a distance.
- On-screen event log with date & time stamp provides trackable infusion history.
- Four infusion modes on AJ5808A numeric keypad version.
- Red light alarm with six different alarm messages. Two different acoustic alarms.



Dual pole clamp enables two infusion pumps on the same level to work as a dual-channel infusion pump; with independent programming on each channel.



• Optional dual pole clamp



### Easy to use, with improved safety features.

AJ5808 and AJ5808A volumetric infusion pumps are designed for a wide range of infusion treatments including general IV infusions, antibiotics, blood infusions, enteral infusions, parenteral nutrition etc. The safety features include anti-free-flow system, on-screen event log with date & time stamp, automatic keypad lock and adjustable occlusion pressure. They can be used in hospital, surgical centers, homecare, or ambulance with the built in rechargeable NiMH battery or plug-in power supply.

Both models can be conveniently connected with our FW-300 fluid warmer for relief of pain and discomfort associated with infusion of cold fluid.



**FW 300  
FLUID WARMER**

